

# FiveGo pH meter F2

Simplify your Mobile Lab. *FiveGo Portable F2 pH/mV Meter*

Measurement made easy

A large, well-structured display, intuitive button arrangement and simple menu ensure measurements can be performed in just a few clicks.

Ergonomic design for one-handed operation

Thanks to its ergonomic design, handling the FiveGo is a comfortable experience. The foldable stand provides optimal viewing on flat surfaces.

IP67 Water and dustproof

The IP67 waterproof protection makes the FiveGo portables able to withstand wet and demanding environments.

Material No.: 30266946



## Specifications - FiveGo pH meter F2

Parameter	pH / mV
pH measuring range	0 – 14
pH resolution	0.01
pH accuracy (±)	0.01
mV measuring range	-1,999 – 1,999
mV resolution	1
mV accuracy (±)	1
Temperature Range	0 °C – 100 °C
Temperature Resolution	0.1 °C
Temperature accuracy (±)	pH: ± 0.5 °C