



VEVOR Digital Rotary Viscometer Viscosity 10-100000 mPa.s Viscosity Meter Tester NDJ-5S Testing Machine Display Viscometer Digital LCD Display

Product features

Digital Viscometer Measurement range: $1-1 \times 100000 \text{ mPa}\cdot\text{s}$; Relation humidity: $\leq 80\%$; Ambient temperature: $5^{\circ}\text{C}-35^{\circ}\text{C}$

Viscosity Meter Tester Rotor specification: 1#,2#,3#,4#,four kind rotors(You can determine viscosity as low as $0.1 \text{ mPa}\cdot\text{s}$ with 0# rotor)

Measurement accuracy: Plus or -2% ; Power: Ac 220 V plus or minus 10% frequency 50 Hz plus or minus 10%

Aluminum material at the bottom of the fuselage, effectively prevent electrostatic interference

NDJ-5S digital rotary viscometer has been widely used to determine and measure the liquid viscosity in many applications such as grease, painting, pharmacy and adhesives