

# Scilogex Levo Plus



## Product Highlights

Scilogex Levo Plus Motorized Pipette Filler, includes 2 x 0.45um Filters, AC Charger & Wall/Bench Stand, 110-220V, 50/60Hz, US plug

Levo Plus is an easy to use pipette controller for use with 1-100ml serological pipettes. Simply attach any size serological pipette into the silicone holder and you are ready to pipette. The LCD display shows the battery level and with adjustable coarse speed setting. Touch sensitive triggers control the fine speed of aspiration and pipetting process. A replaceable 0.45um hydrophobic filter prevents liquids from contaminating the inside of the unit.

CD Display shows battery charge and coarse speed setting

Lightweight ergonomic design, comfortable concave triggers

Powerful lithium-ion battery, 8 hour continuous use per charge

Can fill a 25ml pipette in 5 seconds

Touch sensitive triggers for fine speed control

Gravity dispense, with blow-out for TD pipettes

Autoclavable nose-cone

Can be used while charging, no down-time

Model Name Scilogex Levo Plus : 74020002

SKU : DMPF000

Manufacturer : Scilogex

Manufacturer ID : 74020002

Brand : Scilogex