

DT5-6A low-speed centrifuge



DT5-6A low-speed centrifuge

Features:

Brushless Motor. Digital display. Accurate speed control. Preset rotation speed, time. Safety lid lock, the motor will stop once open the lid. New panel design, easy to operate. A preferable separation device for medical health, laboratory, food, environmental protection, scientific research and teaching.

Technical data	DT5-6A
Max Speed	4000r/min
Max RCF	1780×g
Time Range	0~99min
Power Supply	AC110V - AC220V
Package dimension	340 mm x 375 mm x 270 mm
Weight	13kg
Capacity	10/15ml×8 Angle rotor