



## รายละเอียดสินค้า : รุ่น DR701

1. Range: 0.0 to 50.0% Brix (ATC)
2. TDS : 0.0-25.0% (ATC)
3. Refractive Index : 1.3330-1.4200nD
4. Resolution :  $\pm 0.1\%$  Brix ;  $\pm 0.0001$ nD Refractive Index
5. Accuracy :  $\pm 0.2\%$  Brix ;  $\pm 0.0003$ nD Refractive Index
6. Automatic Temperature Compensation (ATC) : between 10 and 40°C (50 - 104°C)
7. Measurement Time : Approximately 1.5 seconds
8. Minimum Sample Volume : 2 metric drops
9. Sample Cell : Stainless steel ring and flint glass prism
10. Case Material : ABS
11. Enclosure Rating : IP65
12. Battery Type / Life : (1) 9V / 5000 readings
13. Auto-Off : After 3 minutes of non-use
14. Dimensions : 58\*25\*121mm
15. Weight : 350g