

# HEATING MANTLE

---



## Heating Mantle



Intelligent temperature control



Heating evenly, heating up fast

### TECHNICAL DATA

Model	HM		HMC		HMS			HMSC	
Capacity	100ml	250ml	500ml	1L	2L	3L	5L	10L	20L
Function	Heating		Heating		Heating/magnetic stirrer			Heating/magnetic stirrer	
Heating Power	100W	150W	250W	350W	450W	600W	800W	1200W	2400W
Temperature Control	/		Digital temperature control			/		Digital temperature control	
Max Temperature	Room temp - 380°C		Room temp - 300°C		Room temp - 380°C			Room temp - 300°C	