

## Digital Hand-held "PEN" Refractometer PEN-PRO


 Enlargement

 Quantity  

The PEN-PRO has a measurement range of Brix 0.0 to 85.0% which can be used to measure almost any sample. There are 2 ways to measure with the PEN-PRO. Either dip the tip into the sample and press the START key OR press the START key and touch the tip into the sample. There is also the new feature to perform continuous measurements. The temperature range of the "Dip-Type" Refractometer PEN-PRO has been improved from 60°C to 100°C. No more wait ! Measure hot samples while in the cooking process, such as soups and sauces.

#### Specifications

Model	PEN-PRO
Cat.No.	3730
Range	Brix : 0.0 to 85.0%
Accuracy	± 0.2%
Temperature compensation	10 to 100°C
Measurement Time	2 seconds
Power Supply	1 × Size AAA alkaline battery

#### Options

- USB BATTERY CHARGER : RE-34801
- STRAP : RE-39410
- \* Calibration Certificate : Contact an ATAGO representative for details.

Patent for Design Registration No. ZL200830008615.4 (China), D127363 (Taiwan)  
 Patent pending in all other countries around the world.