

## Premium Bunsen Burner with Flame Stabilizer and Gas Adjustment - Liquid Propane (LP)

Model : CH0094LP

- Premium bunsen burners with air flow and gas adjustments
- Overall height is 6.5"
- Accepts 6mm to 8mm inlet tubing

Includes: (1) Bunsen Burner

These premium bunsen burners allow for complete control with both air mix adjustment and gas flow adjustment. The flame stabilizer makes them safer than standard bunsen burners. Accepts 6mm to 8mm inlet tubing.

Burner Specifications: Tube Dia. - 13 mm, Height - 160 mm, Weight - 160 g.

