

## pH Electrode # ECFC7252101B



### Specifications

Product Specification	Description
Parameter	pH
Range	0 to 13 pH
Temperature Range	0 to 80 °C
Liquid Junction Type	Porous HDPE pin
Internal Reference Type	Ag/AgCl
Sealed/Refillable	Sealed
Reference Junction	Single
Refilling Reference Electrolyte	–
Dimensions (Shaft)	90 x 12 mm
Cable Length	1 m
Connector	BNC
Description	General purpose plastic-body pH electrode
Used With	All pH meters with BNC input connector