

Arioso Series

Precision Elegance

Efficiently Innovative System of Human Corporation.


• Extensible point of use



• System Display



POWER (RUN)



CALIBRATION

New model, ARIOSO is the world's top which is embedded with state of the art technology. It will satisfy you beyond your imagination with the world's first developed 2-in-1 system (producing Pure water and Ultra pure water from 1 system), controlled by a large touch screen, meeting the various applications, elegant and practically scientific design.

• Features

- Certified CE, ISO 9001, ISO 14001
- Awarded Korean excellent products
- Auto calibration for water sensor
- Extensible point of use
- Frost alarm message
- 3, 5 steps water level sensor
- Filter exchange indicator
: Beep & character display
- Languages : English, Korean, etc
- Self-diagnosis functions

• Applications (Type I Grade Water)

- Analytical instruments
- AAS, ICP/MS, IC
- HPLC, GC
- TOC
- Cell & tissue culture
- In vitro fertilization
- Electrophoresis
- RNase, DNase & DNA free
- Molecular biology



• Arioso Power Series (RO & UP System / 2-in-1 System)

Type III & I Grade Water

	Arioso Power I II III	Arioso Power +70 I II III	Arioso Power (TOC) II III	Arioso Power +70 (TOC) I II III
Product flow rate (RO, UP)	Max. 15 25 35 L/Hr	Max. RO : 15 25 35 L/Hr UP : 70 L/Hr	Max. 25 35 L/Hr	Max. RO : 15 25 35 L/Hr UP : 70 L/Hr
Display & Controller	Large color 70 x 52 touch screen, Pure Water (RO) : 0.2 to 250.0µs/cm, Ultra Pure (UP) : Up to 18.2MΩ·cm			

• RO System (Pure System)

Type III Grade Water

	Arioso RO 180	Arioso RO 280	Arioso RO 380
Product flow rate (RO)	Max. 15 L/Hr	Max. 25 L/Hr	Max. 35 L/Hr
Display & Controller	Large color 70 x 52 touch screen, Pure Water (RO) : 0.2 to 250.0µs/cm		

• UP System

Type I Grade Water

	Arioso UP 900	Arioso UP 900 (TOC)
Water quality display	Ambient temp : -10 to 100°C, Feed water temp : 0 to 100°C	←
	Feed water : 0 to 990µs/cm, Ultra pure (UP) : Up to 18.2MΩ·cm	←
	-	TOC : 0 to 200 ppb
Product flow rate (L/min)	1.5	1.5
Display & Controller	Large color 70x52 touch screen, , Ultra Pure (UP) : Up to 18.2MΩ·cm	

Optional Accessories	<ul style="list-style-type: none"> • Dispenser • Silver A/C filter • Vent filter, UV lamp for PE water tank 	<ul style="list-style-type: none"> • Volumetric flow meter • Low pressure valve • Level sensor (3, 5 steps) 	<ul style="list-style-type: none"> • Printer • Wall-mounting bracket • Frame 	<ul style="list-style-type: none"> • Foot-Pedal switch • Water tank (20L, 30L, 40L, 60L) 	<ul style="list-style-type: none"> • Pretreatment system • Ion filter • Water leak sensor
----------------------	--	--	---	--	--

Arioso Dispenser Unit

Efficient Dispenser for Water Purification System

Control the dispensing Ultra Pure (UP) water efficiently with Touch Screen

• Features

- Dispenser : Max. 3m
- Dispensing mode
 - Low flow
 - Medium flow
 - High flow
 - Continuous flow
- Error Display
 - Filters Replacement Period
- Volumetric & Timer Dispensing
 - : Self calibration
- Can be connected to Arioso series
- Convenient and optimized dispenser to produce type I grade water.
- Built-in various functions such as dispensing speed (3 steps), volumetric & timer dispensing.
- Display of error, filter replacement period and so on.

• Stand type



Arioso Series




Dispenser Unit

• Dispenser Display



• On the wall type



• Built-in instrument type



• Water tank (Pedestal option)



Human Corporation

Human Bldg, 36, Garak-ro (225-10, Seokchon-dong),
Songpa-gu, Seoul, Korea
Tel. +82-2-418-4182 Fax. +82-2-423-6430
human@humancorp.co.kr www.humancorp.co.kr



Distributor