

Product Name : Mercury (II) Sulfate, Grade AR

Product Code : M2038-1

CAS Number : 7783-35-9

Formula : HgSO_4

M. wt : 296.65 g/mol

Packing : 100 g, 250 g

Parameter	Specification
Assay (complexometric)	min. 99 %
Chloride (Cl)	max. 0.003 %
Nitrates (NO_3)	max. 0.005 %
Cadmium (Cd)	max. 0.0001 %
Calcium (Ca)	max. 0.003 %
Copper (Cu)	max. 0.0005 %
Iron (Fe)	max. 0.001 %
Lead (Pb)	max. 0.0005 %
Magnesium (Mg)	max. 0.003 %
Nickel (Ni)	max. 0.001 %
Potassium (K)	max. 0.002 %
Zinc (Zn)	max. 0.001 %
Sulfated ash after reduction	max. 0.02 %
Mercury (I) (as Hg)	max. 0.05 %
Suitability for COD	passes test

This certificate is an electronic copy and does not require signature